

# Seronorm™ Immunoassay

The all-in-one immunoassay control since 2002



## Analytes

17- $\alpha$ -OH-Progesterone	CK-MB	IGF-1	SHBG
25(OH) Vitamin D	Cortisol	IgFBP-3	T3, free
25(OH) Vitamin D3	C-peptide	Insulin	T3, total
AFP	DHEA-Sulfate	LH	T4, free
Aldosterone	Digoxin	Methylmalonic acid	T4, total
AMH	EPO	Myoglobin	TBG
Androstenedione	Estradiol	NT-proBNP	Testosterone
Anti-Tg	Ferritin	Ostase	Theophylline
Anti-TPO	Folate	Progesterone	Thyroglobulin
$\beta$ 2-microglobulin	FSH	Prolactin	Troponin I
$\beta$ -hCG, total	hCG, total	Prostatic acid phosphatase	Troponin T
CA 125	HE4	PSA, complexed	TSH
CA 15-3	hGH	PSA, free	Vitamin B12
CA 19-9	Homocysteine	PSA, total	
CEA	IgE	PTH, intact	

## Main features & benefits



**The “all-in-one” immunoassay control**  
Limits the number of different controls needed in your lab, saves time, space and money.



**Based on 100% human matrix**  
No stabilizer and no preservative contribute to a high level of patient-likeness.



**Comprehensive product documentation for four different levels**  
Provides analytical data for more than 50 analytes for clinically relevant concentrations.



**Independent QC material**  
Provides an objective assessment of the entire analytical system



**Liquid version**  
No need for reconstitution. Eliminates potential handling errors. Ready to use right after thawing. Can be refrozen.



**Lyophilized version**  
4 years shelf life reduce the frequency of lot change: resource and time-related consequences

## Stability

Liquor	Shelf life		Open vial stability	
	≤ -20 °C	3 years Δ	≤ -20 °C	2-8 °C

  

Lyophilized	Shelf life		Open vial stability	
	2-8 °C	4 years Δ	≤ -20 °C	2-8 °C

- Δ With a few exceptions
- \* When refrozen within 30 minutes
- \*\* When reconstituted and frozen within 30 minutes

## Order Information

Art. no	Product	Size
207005	Seronorm™ Immunoassay Liq L-1*	12 x 3 mL
207105	Seronorm™ Immunoassay Liq L-2*	12 x 3 mL
207205	Seronorm™ Immunoassay Liq L-3*	12 x 3 mL
206005	Seronorm™ Immunoassay Lyo L-1*	12 x 3 mL
206105	Seronorm™ Immunoassay Lyo L-2*	12 x 3 mL
206205	Seronorm™ Immunoassay Lyo L-3*	12 x 3 mL
207305	Seronorm™ Immunoassay Liq. Low*	12 x 3 mL

\* Not available in the US